AxSYM Active-B12

(holotranskobalamina)

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1. JNrn Schneede, Prerequisites for establishing general recommendations for diagnosis and treatment of vitamin B₁₂ deficiency and cost-utility evaluation of these guidelines. *Scand J Clin Lab Invest* 2003;63:369-376.
2. Nexo E, Christensen A-L, Hvas A-M, Petersen TE, Fedosov SN. Quantification of holo-transcobalamin, a marker of vitamin B₁₂ deficiency. *Clin Chem* 2002;48:561-562.
3. Ulleland M, Eilertsen I, Quadros EV, Rothenberg SP, Fedosov SN, Sundrehagen E, Orning L. Direct assay for cobalamin bound to transcobalamin (Holo-transcobalamin) in serum. *Clin Chem* 2002;48:526-532.
4. Lindgren A, Kilander A, Bagge E, Nexo E. Holotranscobalamin – a sensitive marker of cobalamin malabsorption. *Eur J Clin Invest* 1999;29:321-329.
5. Tisman G, Vu T, Amin J, Luszko G, Brenner M, Ramos M *et al.* Measurement of red blood cell vitamin B₁₂: a study of the correlation between intracellular B₁₂ content and concentrations of plasma holotranscobalamin II. *Am J Hematol* 1993;43:226-229.
6. Herbert V, Fong W, Gulle V, Stopler T. Low holotranscobalamin II is the earliest serum marker for subnormal vitamin B₁₂ (cobalamin) absorption in patients with AIDS. *Am J Hematol* 1990;34:132-139.
7. Obeid R *et al.* Cobalamin status (holo-transcobalamin, methylmalonic acid) and folate as determinants of homocysteine concentration. *Clin Chem* 2002;48:2064-2065.
8. Herrmann W *et al.* Functional vitamin B₁₂ deficiency and determination of holotranscobalamin in populations at risk. *Clin Chem Lab Med* 2003;41:1478-1488.
9. Herrmann W *et al.* Vitamin B₁₂ status, particularly holotranscobalamin II and methylmalonic acid concentrations, and hyperhomocysteinemia in vegetarians. *Am J Clin Nutr* 2003;78:131-167.
10. Lloyd-Wright Z *et al.* Holotranscobalamin as an indicator of dietary vitamin B₁₂ deficiency. *Clin Chem* 2003;49:2076-2078.
11. Refsum H *et al.* Hyperhomocysteinemia and elevated methylmalonic acid indicate a high prevalence of cobalamin deficiency in Asian Indians. *Am J Clin Nutr* 2001;74:233-241.
12. Refsum H, Smith AD. Low vitamin B₁₂ status in confirmed Alzheimer's disease as revealed by serum holotranscobalamin. *J Neurol Neurosurg Psychiatry* 2003;74:959-961.
13. Chen X *et al.* Influence of cobalamin deficiency compared with that of cobalamin absorption on serum holotranscobalamin II. *Am J Clin Nutr* 2005;81:110-114.

Abbott GmbH & Co. KG
Abbott Diagnostics Europe
Max-Planck-Ring 2
65205 Wiesbaden
Germany
Tel. (+49) 61 22 58 0
Fax (+49) 61 22 58 12 44
www.abbottdiagnostics.com

Axis-Shield plc
The Technology Park
Dundee DD2 1XA
UK
Tel. (+44) 1382 422 000
Fax (+44) 1382 422 088
www.axis-shield.com